In vitro multiplication of sweetpotato - OP59

TITLE	In vitro multiplication of sweetpotato - OP59
OWNER	Head_IVU,R.Vollmer
APPROVER	Head Genebank
APPROVAL DATE	
ISSUE DATE	Oct 20, 2010 15:08
CONTRIBUTORS	
CITATION	
KEYWORDS	accredited , procedure , head_ivu , head_genebank , r_vollmer , op59
NEXT REVISION	
DOCUMENT ID	OP59
VERSION NUMBER	53

INTERNATIONAL POTATO CENTER - CIP			
			
Version	Date	Author	Comment
53	Oct 20, 2010 15:08	Rainer Vollmer	minor changes - formatting
52	Oct 20, 2010 14:23	Rainer Vollmer	formatting
51	Sep 16, 2010 10:14	Rainer Vollmer	orthographical corrections
50	Sep 14, 2010 08:57	Rainer Vollmer	
49	Sep 09, 2010 14:05	Rainer Vollmer	
48	Sep 08, 2010 15:10	Rainer Vollmer	minor format changes
47	Sep 08, 2010 08:56	Rainer Vollmer	minor changes
46	Sep 08, 2010 08:37	Rainer Vollmer	
45	Sep 07, 2010 19:56	Luis Rojas	Authorized by Rainer Vollmer labels & owner
44	Sep 07, 2010 18:35	Luis Rojas	Authorized by Rainer Vollmer

INTRODUCTION

This technique is based on the principle that appropriate culture conditions induce the growth of pre-existent terminal or axillar buds, resulting in a new plantlet. The nutritional and hormonal conditions of the medium break the bud's latency and promote its rapid growth. Callus formation, followed by shoot regeneration, must be avoided because it can affect the genetic integrity of the clone. Following the protocol described below, after 4-5 weeks of *in vitro* culture plantlets may be conserved, multiplied again, or distributed to the national, international y CIP internal consignees.

SCOPE

This procedure covers the multiplication of *in vitro* sweetpotato germplasm for distribution and conservation purposes. The multiplication procedure is realized by a properly qualified technician (or auxiliary technician).

SAFETY

The operative personnel should know the following protocols:

- Preparation of culture media - OP74 (MPB culture medium)
- Good practices during *in vitro* bank activities

MATERIALS

EQUIPMENT	
Autoclave	Medium dispenser
pHmeter	Analytical balance
Laminar flow chamber	Magnetic Stirrer
Refrigerator	Sterilizer
Environmental monitoring device (HOBO Logger)	Microwave oven
OTHER MATERIALS	
Glass test tubes (13x100 mm, 18x150 mm, 25x150 mm)	Cotton
Test tube caps	Labels
Sterile petri dishes	Wash bottle with alcohol (70°)
Forceps	Burner (with alcohol 96°)
Scalpels	Sterile paper sheets
Blades No.10	Racks
Saran wrap or Parafilm	Plastic Beakers (1, 2 and 4 litres)
Spatula	Magenta vessels
Culture medium (MPB)	Glass vessel with alcohol (96°)

PROCEDURE

1. Clean the flow chamber according to the "[Good practices during in vitro bank activities](#)" protocol. Set up the following materials inside the flow chamber: sterile paper or petri dishes, scalpels with blades, forceps, piece of cotton, alcohol, burner and sterilizer.
2. Take the rack with the accessions to be multiplied to the transference room. Register the following information in the data base: type of medium, batch of medium, purpose of propagation (activity), name of technician, future rack location. Print out the labels using the barcode system according to the number of vessels to be propagated.
3. Under aseptic conditions, the plantlets are taken out from the *in vitro* vessel and placed onto a steril petri dish or paper sheet using sterilized forceps (**Figures 1A and 1B**). It is recommended to use 4-5 week-old plantlets grown in MPB propagation medium (**Table 1**) (see [Preparation of culture media - OP74](#) protocol). Each accession is multiplied separately.
4. Using a sterilized scalpel with a No. 10 blade, remove the tip, leaves and roots (**Figure 1C**) and cut the stem into several segments of approximately 1.0 - 1.5 cm, each containing 2 to 3 buds (**Figure 1D**).
5. Colocate one explant per 13x100 mm tube (national and international distribution), two per 25x150 mm tube (sub-culture, and CIP internal distribution), two per 18x150 mm tube (conservation activity), or five per magenta vessel (CIP internal distribution). Place the explants on the surface of fresh-sterilized MPB propagation medium, maintaining the bud facing upwards (**Figure 1E**).
6. Discard the paper sheet and the plant remains after the multiplication of each accession. Cap and seal the tubes or magenta vessels with a gas-permeable plastic tape (Saran wrap or Parafilm) and label correctly (**Figure 1F**). Use for each accession a new paper sheet. Sterilize the tools after each accession placing them in the sterilizer for at least 30 seconds.
7. Cultures are transferred to the culture incubation room, in which can be found the following controlled physical environment conditions : temperature of 23-25±1 °C, photon flux density (PFD) of 45 mol/(m².seg), with a 16h/8h photoperiod of light and darkness. The environment of the culture incubation chamber is monitored according to a predefined programme. Equipment performance is also checked. The details of environmental and equipment monitoring can be found at the following protocol: [Equipment and environment monitoring \(OP71\)](#). (*)
8. After 5 days of culture, visually evaluate the plantlets for bacteria/fungi contamination, also verify correct labeling, correct number of vessels per accession, and correct number of plants per vessel. Record the evaluation data in the IVU-database. Plantlets showing signs of contamination are discarded and sterilized by autoclaving.
9. After 14 days of culture, visually evaluate the viability of the plantlets (**Table 2**). Record the evaluation data in the IVU-database. Plantlets with abnormal growth are discarded and sterilized by autoclaving.
10. Distribution activity: Perform on delivery date to ADU a visual quality check of the plant material. Record the evaluation data in the IVU-database. Plantlets with normal development are ready for their transfer to ADU.
(*) Illuminance is constantly checked using an environmental monitoring device (HOBO logger). Fluorescents with low illuminance capacity are replaced by new ones (after 3,000 h of use approximately).


Figure 1. Procedure of sweetpotato *in vitro* multiplication. **A)** Extraction of material. **B)** Colocation of the mother plant(s) on petri dish o steril paper. **C)** Removing of tip, leafs, and roots. **D)** Cutting of segments (with 2-3 buds per segment). **E)** Colocation of the explants in fresh MPB medium. **F)** Sealing and labelling of the tube(s)

Table 1. Composition of MPB medium

Ingredient	Quantity (per litre)
MS salts*	4.33 g
Ascorbic acid	200 mg
Calcium nitrate	100 mg
Calcium panthotenate	2 mg
L-Arginine	100 mg
Putrescine-HCl	20 mg
Sucrose	30 g
Phytigel	3 g
pH	5.7
*Murashige and Skoog basal Salts (1962) Supplier: CAISSON	

Table 2. Indicators for viability of *in vitro* sweetpotato plantlets.

General Aspect	Chlorosis	Root formation	Necrosis	Vitrification
Good	Absent	Normal	Absent	Absent
Intermediate	Incipient	Intermediate	Incipient	Intermediate
Regular	Intermediate	Poor	Intermediate	High
Dead	High	Absent	High	

INTERNAL QUALITY CONTROL

1. Quality check list (distribution activity).
2. Table of distribution times (distribution activity).
3. Automatic notification 7 days before sending to ADU (distribution activity).
4. Colocation of original and new label onto one of the multiplied tubes (conservation activity).

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